Manual de instrucciones Bedienungsanleitung

pipet4u® pro

Pipeta de microvolumen

Table of Contents

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General information - The **pipet4u® pro** pipettes are piston-operating pipettes for accurate and precise sampling and dispensing of liquid volumes in the range from 0.1 µL to 10 mL. They operate according to the air displacement principle and are produced and tested according to EN ISO 8655. The **pipet4u® pro** pipettes are available as single channel pipettes with variable and fixed volume and multi channel pipettes with variable volume.


Intended use - The **pipet4u® pro** pipettes are general laboratory pipettes.

Safety - Before using the device, the instruction manual must be read. Follow the safety instruction, wear protective equipment while using the device, especially when working with infectious and dangerous samples. Work in such a way that neither user nor other persons are endangered. Only use the pipette with conform tips. Check the proper condition of the device before each use.

Delivery Package


1 pc. Pipette
1 pc. Calibration / disassembling tool
1 pc. Tool to remove colour-coding drop
1 pc. Silicone grease
1 bag Drop for colour coding
1 bag Sample tips
1 bag O-rings (available for multi channel pipettes)
1 pc. Instruction manual
1 pc. Warranty card
1 pc. Calibration report according to EN ISO 8655

Calibration/disassembling tool - The enclosed tool is designed for easy and quick maintenance and calibration of the pipette.

Pipette identification - Each **pipet4u® pro** pipette can be individually colour coded, using the included tool for drop removing.

Pipette Description

pipet4u® pro

Single channel pipette

1 Push button

Used to aspirate and dispense the liquid and to set the volume (variable volume pipettes)

2 Top closure

Used to adjust the calibration with the help of the calibration tool

3 Drop

The drop is colour-coded as per the volume of the pipette for easy identification

4 Tip ejector

Tip ejector moves the ejector collar and ejects the pipette tip

5 Volume print

Indicates the volume of the pipette

6 Main body

7 Ejector collar

Used to eject the pipette tip

8 Tip cone

The part where the pipette tip is fitted

9 Volume display

Shows the volume set with the push button (variable volume pipettes) and the nominal volume of fixed volume pipette

10 Body Grip

Prevents transfer of body heat to internal components to ensure accurate pipette operation and also provides proper grip while holding the pipette

pipet4u® pro

Multi channel pipette

Upper part of the multi channel pipette corresponds to the single channel pipette as described above.


1 Ejector collar

Used to eject the pipette tips

2 Tip cone


The part where the pipette tips are fitted

Sealing and ejecting tips


Before fitting the tip(s) make sure that the tip cone(s) is/are clean. Use the correct tip according to the volume range or the colour code. Press the tip(s) firmly onto the tip cone(s) to ensure an airtight seal¹. The pipettes are equipped with a tip ejector to help remove the tip(s) and to avoid safety hazards associated with contamination. The tip ejector needs to be pressed downwards to ensure proper tip ejection². (For 10000 μ L and 2000-10000 μ L pipettes tip ejection function is not available).

Volume setting


The volume of the **pipet4u[®]pro** pipette is shown through the hand-grip window¹. For variable pipettes, the desired volume is set by turning the push button clockwise or anticlockwise². Using excessive force to turn the push button outside the range may jam the mechanism and damage the pipette.

Forward pipetting

Press the push button to the first stop **1**. Dip the tip(s) attached to the pipette vertically 2-3 mm into the liquid and smoothly release the push button **2**. Withdraw the tip(s) slowly from the liquid, touching against the wall of the vessel to remove any remaining liquid. Place the tip(s) on the wall of the vessel at an angle. Dispense the liquid into the receiving vessel by gently pressing the push button to the first stop **3**. After a short delay press the push button to the second stop **3**. This blow-out will empty the tip(s) completely and ensure accurate pipetting. Release the push button to the ready position **4**.

Reverse pipetting

Reverse pipetting is recommended for viscous solutions, i.e. solutions which have a tendency to foam, or for dispensing very small volumes. Press the push button to the second stop **1**. Dip the tip(s) attached to the pipette vertically 2-3 mm into the liquid and smoothly release the push button **2**. This will fill the tip(s) with a additional volume larger than the set volume. Withdraw the tip(s) slowly from the liquid, touching against the wall of the vessel to remove any excess liquid. Dispense the liquid into the receiving vessel by gently pressing the push button to the first stop **3**. This volume is equal to the set volume. Hold the push button in this position. The liquid that remains in the tip should not be dispensed. The remaining liquid can be discarded with the tip or delivered back into the vessel with the original solution **4**. Release the push button to the ready position **5**.


Pipetting recommendations - Only aspirate liquid into the pipette when a tip(s) is/are attached. Hold the pipette vertically when aspirating and dip the tip(s) only 2-3 mm into the liquid. Pre-rinse the tip(s) 5 times before aspiration by filling and emptying the tip(s). Always control the push button movements with the thumb for consistency. Allow liquid, tip(s) and pipette to equilibrate to the ambient temperature before use. Do not keep the pipette in your hand while not working to avoid transferring body heat, which may result in incorrect dispensing volume. Use the correct pipette tip(s) designed for the use with the particular pipette. Select the correct pipetting technique (e.g. reverse, forward pipetting) depending on the nature of the liquid. Do not lay the pipette down with a filled tip, otherwise liquid may flow into the pipette and contaminate it.

Maintenance- It is recommended to check the performance of your pipette regularly e.g. every 3 months and always after in-house service or maintenance. To maintain the best results from your pipette, each unit should be checked every day for cleanliness. Pay special attention to the tip cone(s).

Disassembling/reassembling and cleaning

To clean and decontaminate or replace the lower parts of the single channel pipette, proceed as follows:

- 1 Hold down the tip ejector.
- 2 Place the tooth of the opening tool between tip ejector and tip ejector collar to release the locking mechanism.
- 3 Carefully release the tip ejector and remove the ejector collar.
- 4 Place the wrench end of the opening tool over the tip cone and turn it clockwise.
- 5 After removing the tip cone, wipe the piston, the a-ring and the tip cone with ethanol and a lint free cloth. (cleaning step)
- 6 Before replacing the tip cone grease the piston slightly using the silicon grease provided. (cleaning step)
- 7 After reassembling use the pipette several times to make sure that the grease is spread evenly.
- 8 Check the calibration of the pipette.

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

Top assembly parts - pipet4u

Pipetas monocanal y multicanal


Description	
1 Push Button	8 Washer
2 Top closure	9 Shaft
3 Spring stopper	10 Bottom closure snap fitting
4 Primary spring support	11 Tip ejector
5 Spring primary	12 Spring tip ejector
6 Main body	13 Drop
7 Counter	14 Tip ejector pusher support

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
Pipetas multicanal

Description

- | | |
|---|------------------------|
| 1 | Piston housing |
| 2 | Steel piston |
| 3 | Tip cone housing |
| 4 | Sealing O-rings |
| 5 | Tip cone outer O-rings |

Instruction Manual

Specifications

Single channel pipettes - variable volume

Volume Range	Increments	Test Volume	Inaccuracy (±)%	Imprecision (±)%
0.1 µL - 2.5 µL	0.01 µL	0.25 µL	12.00%	6.00%
		1.25 µL	3.00%	3.00%
		2.5 µL	2.50%	1.60%
0.5 µL - 10 µL	0.1 µL	1 µL	2.50%	1.50%
		5 µL	2.00%	1.00%
		10 µL	1.00%	0.80%
2 µL - 20 µL	0.1 µL	2 µL	3.00%	2.00%
		10 µL	1.50%	1.00%
		20 µL	0.90%	0.40%
5 µL - 50 µL	0.5 µL	5 µL	2.00%	2.00%
		25 µL	0.80%	0.40%
		50 µL	0.60%	0.30%
10 µL - 100 µL	0.5 µL	10 µL	3.00%	1.50%
		50 µL	1.00%	0.50%
		100 µL	0.80%	0.15%
20 µL - 200 µL	1 µL	20 µL	2.00%	0.80%
		100 µL	0.70%	0.30%
		200 µL	0.60%	0.15%
50 µL - 200 µL	1 µL	50 µL	1.00%	0.40%
		100 µL	0.80%	0.20%
		200 µL	0.60%	0.15%
100 µL - 1000 µL	5 µL	100 µL	2.00%	0.70%
		500 µL	1.00%	0.40%
		1000 µL	0.60%	0.20%
200 µL - 1000 µL	5 µL	200 µL	0.90%	0.30%
		500 µL	0.75%	0.25%
		1000 µL	0.60%	0.20%
500 µL - 5000 µL	50 µL	500 µL	2.00%	0.60%
		2500 µL	0.60%	0.30%
		5000 µL	0.50%	0.15%
2000 µL - 10000 µL	100 µL	2000 µL	3.00%	0.60%
		5000 µL	1.20%	0.30%
		10000 µL	0.60%	0.20%

Single channel pipettes - fixed volume

Volume Range	Increments	Test Volume	Inaccuracy (±)%	Imprecision (±)%
5 µL	-	5 µL	1.30%	1.20%
10 µL	-	10 µL	0.80%	0.80%
20 µL	-	20 µL	0.60%	0.50%
25 µL	-	25 µL	0.50%	0.30%
50 µL	-	50 µL	0.50%	0.30%
100 µL	-	100 µL	0.50%	0.30%
200 µL	-	200 µL	0.40%	0.20%
250 µL	-	250 µL	0.40%	0.20%
500 µL	-	500 µL	0.30%	0.20%
1000 µL	-	1000 µL	0.30%	0.20%
2000 µL	-	2000 µL	0.30%	0.15%
2500 µL	-	2500 µL	0.30%	0.15%
5000 µL	-	5000 µL	0.30%	0.15%
10000 µL	-	10000 µL	0.60%	0.20%

8- and 12- channel pipettes - variable volume

Volume Range	Increments	Test Volume	Inaccuracy (±)%	Imprecision (±)%
0.5 µL - 10 µL	0.1 µL	1 µL	4.00%	4.00%
		5 µL	2.50%	2.50%
		10 µL	1.50%	1.50%
5 µL - 50 µL	0.5 µL	5 µL	3.00%	2.00%
		25 µL	1.50%	1.00%
		50 µL	1.00%	0.70%
10 µL - 100 µL	1 µL	10 µL	1.50%	0.75%
		50 µL	1.00%	0.50%
		100 µL	1.00%	0.25%
30 µL - 300 µL	1 µL	30 µL	1.50%	0.75%
		150 µL	1.00%	0.50%
		300 µL	0.80%	0.25%

ahn

AHN Biotechnologie GmbH
Uthleber Weg 14
99734 Nordhausen
Germany

P: +49(0)3631/65242-0

F: +49(0)3631/65242-90

E: info@cappahn.com

www.ahn-bio.de