

ahn maxxLine

Plastic Cryo Boxes

NEW
PRODUCT

Plastic Cryo Boxes with 10x10 grid 130 x 130 mm

10x10 grid – offers 23% more storage and is ideal for laboratories storing large numbers of samples

Durability - high-quality polypropylene (PP)

Temperature resistant - withstand temperatures -90° to +121°C for reliable cryo vial storage and easy autoclaving

Star-foot lock at the bottom of each box well - enabling one-handed tube opening (only with suitable tube) without taking the tube out of the box, fits perfectly with ahn maxxLine cryo vials with star-shaped foot base

Easy sample organisation - moulded alpha-numeric coding on the lid and bottom of each box enable sample location with the box opened or closed


Slanted edges – for easier identification and handling

Bottom vent holes – rapidly increases freezing time and allow liquid drainage

Transparent lids - make it easy to see contents

Height options – 32 mm, 40 mm, 50 mm, 80 mm and 90 mm

Easy sample sorting – with coloured bottoms in blue, green, red, yellow as well as in clear



Plastic Cryo Boxes with 9x9 grid 130 x 130 mm

9x9 grid – fit for most storage applications

Durability - high-quality polypropylene (PP)

Temperature resistant - withstand -90° to +121°C for reliable cryo vial storage and easy autoclaving

Easy sample organisation - moulded alpha-numeric coding on box lid and bottom enable sample location with the box opened or closed

Bottom vent holes – rapidly increases freezing time and allow liquid drainage

Transparent lids - make it easy to see contents

Height options – 50 mm, 80 mm and 90 mm

Easy sample sorting – with coloured bottoms in blue, green, red, yellow as well as in clear

Ordering information

Plastic Cryo Boxes 130x130 mm with 10x10 Dividers

Description	Grid	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine plastic cryo box for 0.5 mL cryo tubes	10x10	130x130x32 mm	blue	case / 36 pcs.	MLX05B
maxxLine plastic cryo box for 0.5 mL cryo tubes	10x10	130x130x32 mm	green	case / 36 pcs.	MLX05G
maxxLine plastic cryo box for 0.5 mL cryo tubes	10x10	130x130x32 mm	clear	case / 36 pcs.	MLX05N
maxxLine plastic cryo box for 0.5 mL cryo tubes	10x10	130x130x32 mm	red	case / 36 pcs.	MLX05R
maxxLine plastic cryo box for 0.5 mL cryo tubes	10x10	130x130x32 mm	yellow	case / 36 pcs.	MLX05Y
maxxLine plastic cryo box for 1.0 mL cryo tubes	10x10	130x130x40 mm	blue	case / 36 pcs.	MLX1B
maxxLine plastic cryo box for 1.0 mL cryo tubes	10x10	130x130x40 mm	green	case / 36 pcs.	MLX1G
maxxLine plastic cryo box for 1.0 mL cryo tubes	10x10	130x130x40 mm	red	case / 36 pcs.	MLX1R
maxxLine plastic cryo box for 1.0 mL cryo tubes	10x10	130x130x50 mm	yellow	case / 36 pcs.	MLX1Y
maxxLine plastic cryo box for 4.0 mL cryo tubes	10x10	130x130x80 mm	blue	case / 36 pcs.	MLX4B
maxxLine plastic cryo box for 4.0 mL cryo tubes	10x10	130x130x80 mm	green	case / 36 pcs.	MLX4G
maxxLine plastic cryo box for 4.0 mL cryo tubes	10x10	130x130x80 mm	clear	case / 36 pcs.	MLX4N
maxxLine plastic cryo box for 4.0 mL cryo tubes	10x10	130x130x80 mm	red	case / 36 pcs.	MLX4R
maxxLine plastic cryo box for 4.0 mL cryo tubes	10x10	130x130x90 mm	yellow	case / 36 pcs.	MLX4Y
maxxLine plastic cryo box for 5.0 mL cryo tubes	10x10	130x130x90 mm	yellow	case / 36 pcs.	MLX5Y
maxxLine plastic cryo box for 1.0 mL cryo tubes	10x10	130x130x40 mm	clear	case / 36 pcs.	MLX1N
maxxLine plastic cryo box for 2.0 mL cryo tubes	10x10	130x130x50 mm	blue	case / 36 pcs.	MLX2B
maxxLine plastic cryo box for 2.0 mL cryo tubes	10x10	130x130x50 mm	green	case / 36 pcs.	MLX2G
maxxLine plastic cryo box for 2.0 mL cryo tubes	10x10	130x130x50 mm	clear	case / 36 pcs.	MLX2N
maxxLine plastic cryo box for 2.0 mL cryo tubes	10x10	130x130x50 mm	red	case / 36 pcs.	MLX2R
maxxLine plastic cryo box, for 2.0 mL cryo tubes	10x10	130x130x80 mm	yellow	case / 36 pcs.	MLX2Y
maxxLine plastic cryo box for 5.0 mL cryo tubes	10x10	130x130x90 mm	blue	case / 36 pcs.	MLX5B
maxxLine plastic cryo box for 5.0 mL cryo tubes	10x10	130x130x90 mm	green	case / 36 pcs.	MLX5G
maxxLine plastic cryo box for 5.0 mL cryo tubes	10x10	130x130x90 mm	clear	case / 36 pcs.	MLX5N
maxxLine plastic cryo box for 5.0 mL cryo tubes	10x10	130x130x90 mm	red	case / 36 pcs.	MLX5R



+49(0)3631/65242-0



contact@ahn-bio.com




www.ahn-bio.com



Plastic Cryo Boxes 130x130 mm with 9x9 Dividers

Description	Grid	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	blue	case / 36 pcs.	CM81-B
maxxLine plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	blue	case / 36 pcs.	CM81-B4
maxxLine plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	blue	case / 36 pcs.	CM81-B5
maxxLine plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	green	case / 36 pcs.	CM81-G
maxxLine plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	green	case / 36 pcs.	CM81-G4
maxxLine plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	green	case / 36 pcs.	CM81-G5
maxxLine plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	clear	case / 36 pcs.	CM81-N
maxxLine plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	clear	case / 36 pcs.	CM81-N4
maxxLine plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	clear	case / 36 pcs.	CM81-N5
maxxLine plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	red	case / 36 pcs.	CM81-R
maxxLine plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	red	case / 36 pcs.	CM81-R4
maxxLine plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	red	case / 36 pcs.	CM81-R5
maxxLine plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	yellow	case / 36 pcs.	CM81-Y
maxxLine plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	yellow	case / 36 pcs.	CM81-Y4
maxxLine plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	yellow	case / 36 pcs.	CM81-Y5



ahn maxxLine

Cardboard Cryo Boxes

Durability - high quality cardboard boxes feature a very hard-wearing, water-repellent coating which protects them from moisture damage and makes them extra resistant to alcohols and mild organic solvents

Temperature resistant - withstand temperatures as low as -80°C without any loss of integrity

Eco-friendly - fully recyclable, environmentally-friendly materials with low carbon footprint


Available to two footprints - 133x133 mm is particularly suited to space-constrained laboratories / 136x136 mm perfect for large sample storage demands

Height options - 32 mm, 50 mm, 75 mm, 100 mm and 130 mm suitable for different tube heights

Easy sample sorting - multiple colours boxes available: white, blue, green, red and yellow

Labeling option - Write-on surface

Optional cardboard dividers - can be used to store a wide range of laboratory items



Ordering information

Cardboard Cryo Boxes 133x133m with dividers

Description	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine cardboard cryo box B9, 32 mm, incl. 9x9 cells	133x133x32 mm	white	case / 36 pcs.	ML22225
maxxLine cardboard cryo box B10, 32 mm, incl. 10x10 cells	133x133x32 mm	white	case / 36 pcs.	ML22226
maxxLine cardboard cryo box B9, 50 mm, incl. 9x9 cells	133x133x50 mm	white	case / 36 pcs.	ML22230
maxxLine cardboard cryo box B10, 50 mm, incl. 10x10 cells	133x133x50 mm	white	case / 36 pcs.	ML22232
maxxLine cardboard cryo box B9, 75 mm, incl. 9x9 cells	133x133x75 mm	white	case / 36 pcs.	ML22235
maxxLine cardboard cryo box B9, 100 mm, incl. 9x9 cells	133x133x100 mm	white	case / 36 pcs.	ML22240
maxxLine cardboard cryo box B9, 130 mm, incl. 9x9 cells	133x133x130 mm	white	case / 36 pcs.	ML22245



+49(0)3631/65242-0



contact@ahn-bio.com




www.ahn-bio.com



Cardboard Cryo Boxes 133x133m without dividers

Description	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine cardboard cryo box, 32 mm	133x133x32 mm	white	case / 36 pcs.	ML22000
maxxLine cardboard cryo box, 32 mm	133x133x32 mm	blue	case / 36 pcs.	ML22001
maxxLine cardboard cryo box, 32 mm	133x133x32 mm	red	case / 36 pcs.	ML22002
maxxLine cardboard cryo box, 32 mm	133x133x32 mm	green	case / 36 pcs.	ML22003
maxxLine cardboard cryo box, 32 mm	133x133x32 mm	yellow	case / 36 pcs.	ML22004
maxxLine cardboard cryo box, 50 mm	133x133x50 mm	white	case / 36 pcs.	ML22005
maxxLine cardboard cryo box, 50 mm	133x133x50 mm	blue	case / 36 pcs.	ML22010
maxxLine cardboard cryo box, 50 mm	133x133x50 mm	red	case / 36 pcs..	ML22015
maxxLine cardboard cryo box, 50 mm	133x133x50 mm	green	case / 36 pcs..	ML22020
maxxLine cardboard cryo box, 50 mm	133x133x50 mm	yellow	case / 36 pcs..	ML22025
maxxLine cardboard cryo box, 75 mm	133x133x75 mm	white	case / 36 pcs..	ML22050
maxxLine cardboard cryo box, 75 mm	133x133x75 mm	blue	case / 36 pcs..	ML22055
maxxLine cardboard cryo box, 75 mm	133x133x75 mm	red	case / 36 pcs..	ML22060
maxxLine cardboard cryo box, 75 mm	133x133x75 mm	green	case / 36 pcs..	ML22065
maxxLine cardboard cryo box, 75 mm	133x133x75 mm	yellow	case / 36 pcs..	ML22070
maxxLine cardboard cryo box, 100 mm	133x133x100 mm	white	case / 36 pcs..	ML22100
maxxLine cardboard cryo box, 100 mm	133x133x100 mm	blue	case / 36 pcs..	ML22101
maxxLine cardboard cryo box, 100 mm	133x133x100 mm	red	case / 36 pcs..	ML22102
maxxLine cardboard cryo box, 100 mm	133x133x100 mm	green	case / 36 pcs..	ML22103
maxxLine cardboard cryo box, 100 mm	133x133x100 mm	yellow	case / 36 pcs..	ML22104
maxxLine cardboard cryo box, 130 mm	133x133x130 mm	white	case / 36 pcs..	ML22080
maxxLine cardboard cryo box, 130 mm	133x133x130 mm	blue	case / 36 pcs..	ML22082
maxxLine cardboard cryo box, 130 mm	133x133x130 mm	red	case / 36 pcs..	ML22084
maxxLine cardboard cryo box, 130 mm	133x133x130 mm	green	case / 36 pcs..	ML22086
maxxLine cardboard cryo box, 130 mm	133x133x130 mm	yellow	case / 36 pcs..	ML22088



Dividers for Cardboard Cryo Boxes (133x133 mm)

Description	Divider Matrix	Divider Height	Sales Unit	Cat. No.
Divider for 32 mm cryo box	7x7 (vial < 17 mm O.D.)	22 mm	case / 36 pcs..	ML22090
Divider for 32 mm cryo box	8x8 (vial < 15 mm O.D.)	22 mm	case / 36 pcs..	ML22091
Divider for 32 mm cryo box	9x9 (vial < 13 mm O.D.)	22 mm	case / 36 pcs..	ML22092
Divider for 32 mm cryo box	10x10 (vial < 12 mm O.D.)	22 mm	case / 36 pcs..	ML22093
Divider for 32 mm cryo box	12x12 (vial < 10 mm O.D.)	22 mm	case / 36 pcs..	ML22094
Divider for 32 mm cryo box	13x13 (vial < 9 mm O.D.)	22 mm	case / 36 pcs..	ML22095
Divider for 32 mm cryo box	14x14 (vial < 8 mm O.D.)	22 mm	case / 36 pcs..	ML22096
Divider for 32 mm cryo box	16x16 (vial < 7 mm O.D.)	22 mm	case / 36 pcs..	ML22097
Divider for 50 mm cryo box	3x3 (vial < 41 mm O.D.)	30 mm	case / 36 pcs..	ML22111
Divider for 50 mm cryo box	4x4 (vial < 31 mm O.D.)	30 mm	case / 36 pcs..	ML22112
Divider for 50 mm cryo box	5x5 (vial < 24 mm O.D.)	30 mm	case / 36 pcs..	ML22113
Divider for 50 mm cryo box	6x6 (vial < 20 mm O.D.)	30 mm	case / 36 pcs..	ML22114
Divider for 50 mm cryo box	7x7 (vial < 17 mm O.D.)	30 mm	case / 36 pcs..	ML22106
Divider for 50 mm cryo box	8x8 (vial < 15 mm O.D.)	30 mm	case / 36 pcs..	ML22108
Divider for 50 mm cryo box	9x9 (vial < 13 mm O.D.)	30 mm	case / 36 pcs..	ML22110
Divider for 50 mm cryo box	10x10 (vial < 12 mm O.D.)	30 mm	case / 36 pcs..	ML22120
Divider for 50 mm cryo box	12x12 (vial < 10 mm O.D.)	30 mm	case / 36 pcs..	ML22115
Divider for 50 mm cryo box	13x13 (vial < 9 mm O.D.)	30 mm	case / 36 pcs..	ML22116
Divides for 50 mm cryo box	14x14 (vial < 8 mm O.D.)	30 mm	case / 36 pcs..	ML22117
Divider for 50 mm cryo box	16x16 (vial < 7 mm O.D.)	30 mm	case / 36 pcs..	ML22118
Divider for 75 mm cryo box	3x3 (vial < 41 mm O.D.)	40 mm	case / 36 pcs..	ML22124
Divider for 75 mm cryo box	4x4 (vial < 31 mm O.D.)	40 mm	case / 36 pcs..	ML22125
Divider for 75 mm cryo box	5x5 (vial < 24 mm O.D.)	40 mm	case / 36 pcs..	ML22126
Divider for 75 mm cryo box	6x6 (vial < 20 mm O.D.)	40 mm	case / 36 pcs..	ML22127
Divider for 75 mm cryo box	7x7 (vial < 17 mm O.D.)	40 mm	case / 36 pcs..	ML22130
Divider for 75 mm cryo box	8x8 (vial < 15 mm O.D.)	40 mm	case / 36 pcs..	ML22132
Divider for 75 mm cryo box	9x9 (vial < 13 mm O.D.)	40 mm	case / 36 pcs..	ML22134
Divider for 75 mm cryo box	10x10 (vial < 12 mm O.D.)	40 mm	case / 36 pcs..	ML22136
Divider for 75 mm cryo box	12x12 (vial < 10 mm O.D.)	40 mm	case / 36 pcs..	ML22137



Dividers for Cardboard Cryo Boxes (133x133 mm)

Description	Divider Matrix	Divider Height	Sales Unit	Cat. No.
Divider for 100/130 mm cryo box	3x3 (vial < 41 mm O.D.)	65 mm	case / 36 pcs..	ML22147
Divider for 100/130 mm cryo box	4x4 (vial < 31 mm O.D.)	65 mm	case / 36 pcs..	ML22148
Divider for 100/130 mm cryo box	5x5 (vial < 24 mm O.D.)	65 mm	case / 36 pcs..	ML22149
Divider for 100/130 mm cryo box	6x6 (vial < 20 mm O.D.)	65 mm	case / 36 pcs..	ML22150
Divider for 100/130 mm cryo box	7x7 (vial < 17 mm O.D.)	65 mm	case / 36 pcs..	ML22140
Divider for 100/130 mm cryo box	8x8 (vial < 15 mm O.D.)	65 mm	case / 36 pcs..	ML22142
Divider for 100/130 mm cryo box	9x9 (vial < 13 mm O.D.)	65 mm	case / 36 pcs..	ML22144
Divider for 100/130 mm cryo box	10x10 (vial < 12 mm O.D.)	65 mm	case / 36 pcs..	ML22146

Cardboard Cryo Boxes 136x136 mm with Dividers

Description	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine cardboard cryo box B9 incl. 9x9 cells	136x136x50 mm	white	case / 36 pcs..	ML22231
maxxLine cardboard cryo box B10 incl. 10x10 cells	136x136x50 mm	white	case / 36 pcs..	ML22233
maxxLine cardboard cryo box w/print inc. 9x9 dividers	136x136x50 mm	white	case / 36 pcs..	ML22933

