

## Plastic Cryo Boxes

Plastic Cryo Boxes with 9x9 grid  
130 x 130 mm

56

**9x9 grid** – fit for most storage applications

**Durability** - high-quality polypropylene (PP)

**Temperature resistant** - withstand -90° to +121°C for reliable cryo vial storage and easy autoclaving



**Easy sample organisation** - moulded alpha-numeric coding on box lid and bottom enable sample location with the box opened or closed

**Bottom vent holes** – rapidly increases freezing time and allow liquid drainage

**Transparent lids** - make it easy to see contents

**Height options** – 50 mm, 80 mm and 90 mm

**Easy sample sorting** – with coloured bottoms in blue, green, red, yellow as well as in clear



# Ordering information

## Plastic Cryo Boxes 130x130 mm with 9x9 Dividers

Description	Grid	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine® plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	blue	case / 36 pcs.	CM81-B
maxxLine® plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	green	case / 36 pcs.	CM81-G
maxxLine® plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	clear	case / 36 pcs.	CM81-N
maxxLine® plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	red	case / 36 pcs.	CM81-R
maxxLine® plastic cryo box for 1.0-2.0 mL cryo tubes	9x9	130x130x50 mm	yellow	case / 36 pcs.	CM81-Y
maxxLine® plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	blue	case / 36 pcs.	CM81-B4
maxxLine® plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	green	case / 36 pcs.	CM81-G4
maxxLine® plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	clear	case / 36 pcs.	CM81-N4
maxxLine® plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	red	case / 36 pcs.	CM81-R4
maxxLine® plastic cryo box for 4.0 mL cryo tubes	9x9	130x130x80 mm	yellow	case / 36 pcs.	CM81-Y4
maxxLine® plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	blue	case / 36 pcs.	CM81-B5
maxxLine® plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	green	case / 36 pcs.	CM81-G5
maxxLine® plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	clear	case / 36 pcs.	CM81-N5
maxxLine® plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	red	case / 36 pcs.	CM81-R5
maxxLine® plastic cryo box for 5.0 mL cryo tubes	9x9	130x130x90 mm	yellow	case / 36 pcs.	CM81-Y5

ahn®maxxLine®

### Cryo Vials

Pair your samples with AHN maxxLine® Cryo Vials and Cryoboxes — designed to work together, built to last.

Please refer page no: 54



## Cardboard Cryo Boxes

58

**Durability** - high quality cardboard boxes feature a very hard-wearing, water-repellent coating which protects them from moisture damage and makes them extra resistant to alcohols and mild organic solvents

**Temperature resistant** - withstand temperatures as low as -80°C without any loss of integrity

**Eco-friendly** - fully recyclable, environmentally-friendly materials with low carbon footprint


**Available to two footprints** - 133x133 mm is particularly suited to space-constrained laboratories / 136x136 mm perfect for large sample storage demands

**Height options** - 32 mm, 50 mm, 75 mm, 100 mm and 130 mm suitable for different tube heights

**Easy sample sorting** - multiple colours boxes available: white, blue, green, red and yellow

**Labeling option** - Write-on surface

**Optional cardboard dividers** - can be used to store a wide range of laboratory items



## Ordering information


### Cardboard Cryo Boxes 133x133m with dividers

Description	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine <sup>®</sup> cardboard cryo box B9, 32 mm, incl. 9x9 cells	133x133x32 mm	white	case / 36 pcs.	ML22225
maxxLine <sup>®</sup> cardboard cryo box B10, 32 mm, incl. 10x10 cells	133x133x32 mm	white	case / 36 pcs.	ML22226
maxxLine <sup>®</sup> cardboard cryo box B9, 50 mm, incl. 9x9 cells	133x133x50 mm	white	case / 36 pcs.	ML22230
maxxLine <sup>®</sup> cardboard cryo box B10, 50 mm, incl. 10x10 cells	133x133x50 mm	white	case / 36 pcs.	ML22232
maxxLine <sup>®</sup> cardboard cryo box B9, 75 mm, incl. 9x9 cells	133x133x75 mm	white	case / 36 pcs.	ML22235
maxxLine <sup>®</sup> cardboard cryo box B9, 100 mm, incl. 9x9 cells	133x133x100 mm	white	case / 36 pcs.	ML22240
maxxLine <sup>®</sup> cardboard cryo box B9, 130 mm, incl. 9x9 cells	133x133x130 mm	white	case / 36 pcs.	ML22245



## Cardboard Cryo Boxes 133x133m without dividers

Description	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine® cardboard cryo box, 32 mm	133x133x32 mm	white	case / 36 pcs.	ML22000
maxxLine® cardboard cryo box, 32 mm	133x133x32 mm	blue	case / 36 pcs.	ML22001
maxxLine® cardboard cryo box, 32 mm	133x133x32 mm	red	case / 36 pcs.	ML22002
maxxLine® cardboard cryo box, 32 mm	133x133x32 mm	green	case / 36 pcs.	ML22003
maxxLine® cardboard cryo box, 32 mm	133x133x32 mm	yellow	case / 36 pcs.	ML22004
maxxLine® cardboard cryo box, 50 mm	133x133x50 mm	white	case / 36 pcs.	ML22005
maxxLine® cardboard cryo box, 50 mm	133x133x50 mm	blue	case / 36 pcs.	ML22010
maxxLine® cardboard cryo box, 50 mm	133x133x50 mm	red	case / 36 pcs..	ML22015
maxxLine® cardboard cryo box, 50 mm	133x133x50 mm	green	case / 36 pcs..	ML22020
maxxLine® cardboard cryo box, 50 mm	133x133x50 mm	yellow	case / 36 pcs..	ML22025
maxxLine® cardboard cryo box, 75 mm	133x133x75 mm	white	case / 36 pcs..	ML22050
maxxLine® cardboard cryo box, 75 mm	133x133x75 mm	blue	case / 36 pcs..	ML22055
maxxLine® cardboard cryo box, 75 mm	133x133x75 mm	red	case / 36 pcs..	ML22060
maxxLine® cardboard cryo box, 75 mm	133x133x75 mm	green	case / 36 pcs..	ML22065
maxxLine® cardboard cryo box, 75 mm	133x133x75 mm	yellow	case / 36 pcs..	ML22070
maxxLine® cardboard cryo box, 100 mm	133x133x100 mm	white	case / 36 pcs..	ML22100
maxxLine® cardboard cryo box, 100 mm	133x133x100 mm	blue	case / 36 pcs..	ML22101
maxxLine® cardboard cryo box, 100 mm	133x133x100 mm	red	case / 36 pcs..	ML22102
maxxLine® cardboard cryo box, 100 mm	133x133x100 mm	green	case / 36 pcs..	ML22103
maxxLine® cardboard cryo box, 100 mm	133x133x100 mm	yellow	case / 36 pcs..	ML22104
maxxLine® cardboard cryo box, 130 mm	133x133x130 mm	white	case / 36 pcs..	ML22080
maxxLine® cardboard cryo box, 130 mm	133x133x130 mm	blue	case / 36 pcs..	ML22082
maxxLine® cardboard cryo box, 130 mm	133x133x130 mm	red	case / 36 pcs..	ML22084
maxxLine® cardboard cryo box, 130 mm	133x133x130 mm	green	case / 36 pcs..	ML22086
maxxLine® cardboard cryo box, 130 mm	133x133x130 mm	yellow	case / 36 pcs..	ML22088



## Dividers for Cardboard Cryo Boxes (133x133 mm)

Description	Divider Matrix	Divider Height	Sales Unit	Cat. No.
Divider for 32 mm cryo box	7x7 (vial < 17 mm O.D.)	22 mm	case / 36 pcs..	ML22090
Divider for 32 mm cryo box	8x8 (vial < 15 mm O.D.)	22 mm	case / 36 pcs..	ML22091
Divider for 32 mm cryo box	9x9 (vial < 13 mm O.D.)	22 mm	case / 36 pcs..	ML22092
Divider for 32 mm cryo box	10x10 (vial < 12 mm O.D.)	22 mm	case / 36 pcs..	ML22093
Divider for 32 mm cryo box	12x12 (vial < 10 mm O.D.)	22 mm	case / 36 pcs..	ML22094
Divider for 32 mm cryo box	13x13 (vial < 9 mm O.D.)	22 mm	case / 36 pcs..	ML22095
Divider for 32 mm cryo box	14x14 (vial < 8 mm O.D.)	22 mm	case / 36 pcs..	ML22096
Divider for 32 mm cryo box	16x16 (vial < 7 mm O.D.)	22 mm	case / 36 pcs..	ML22097
Divider for 50 mm cryo box	3x3 (vial < 41 mm O.D.)	30 mm	case / 36 pcs..	ML22111
Divider for 50 mm cryo box	4x4 (vial < 31 mm O.D.)	30 mm	case / 36 pcs..	ML22112
Divider for 50 mm cryo box	5x5 (vial < 24 mm O.D.)	30 mm	case / 36 pcs..	ML22113
Divider for 50 mm cryo box	6x6 (vial < 20 mm O.D.)	30 mm	case / 36 pcs..	ML22114
Divider for 50 mm cryo box	7x7 (vial < 17 mm O.D.)	30 mm	case / 36 pcs..	ML22106
Divider for 50 mm cryo box	8x8 (vial < 15 mm O.D.)	30 mm	case / 36 pcs..	ML22108
Divider for 50 mm cryo box	9x9 (vial < 13 mm O.D.)	30 mm	case / 36 pcs..	ML22110
Divider for 50 mm cryo box	10x10 (vial < 12 mm O.D.)	30 mm	case / 36 pcs..	ML22120
Divider for 50 mm cryo box	12x12 (vial < 10 mm O.D.)	30 mm	case / 36 pcs..	ML22115
Divider for 50 mm cryo box	13x13 (vial < 9 mm O.D.)	30 mm	case / 36 pcs..	ML22116
Divides for 50 mm cryo box	14x14 (vial < 8 mm O.D.)	30 mm	case / 36 pcs..	ML22117
Divider for 50 mm cryo box	16x16 (vial < 7 mm O.D.)	30 mm	case / 36 pcs..	ML22118
Divider for 75 mm cryo box	3x3 (vial < 41 mm O.D.)	40 mm	case / 36 pcs..	ML22124
Divider for 75 mm cryo box	4x4 (vial < 31 mm O.D.)	40 mm	case / 36 pcs..	ML22125
Divider for 75 mm cryo box	5x5 (vial < 24 mm O.D.)	40 mm	case / 36 pcs..	ML22126
Divider for 75 mm cryo box	6x6 (vial < 20 mm O.D.)	40 mm	case / 36 pcs..	ML22127
Divider for 75 mm cryo box	7x7 (vial < 17 mm O.D.)	40 mm	case / 36 pcs..	ML22130
Divider for 75 mm cryo box	8x8 (vial < 15 mm O.D.)	40 mm	case / 36 pcs..	ML22132
Divider for 75 mm cryo box	9x9 (vial < 13 mm O.D.)	40 mm	case / 36 pcs..	ML22134
Divider for 75 mm cryo box	10x10 (vial < 12 mm O.D.)	40 mm	case / 36 pcs..	ML22136
Divider for 75 mm cryo box	12x12 (vial < 10 mm O.D.)	40 mm	case / 36 pcs..	ML22137



### Dividers for Cardboard Cryo Boxes (133x133 mm)

Description	Divider Matrix	Divider Height	Sales Unit	Cat. No.
Divider for 100/130 mm cryo box	3x3 (vial < 41 mm O.D.)	65 mm	case / 36 pcs..	ML22147
Divider for 100/130 mm cryo box	4x4 (vial < 31 mm O.D.)	65 mm	case / 36 pcs..	ML22148
Divider for 100/130 mm cryo box	5x5 (vial < 24 mm O.D.)	65 mm	case / 36 pcs..	ML22149
Divider for 100/130 mm cryo box	6x6 (vial < 20 mm O.D.)	65 mm	case / 36 pcs..	ML22150
Divider for 100/130 mm cryo box	7x7 (vial < 17 mm O.D.)	65 mm	case / 36 pcs..	ML22140
Divider for 100/130 mm cryo box	8x8 (vial < 15 mm O.D.)	65 mm	case / 36 pcs..	ML22142
Divider for 100/130 mm cryo box	9x9 (vial < 13 mm O.D.)	65 mm	case / 36 pcs..	ML22144
Divider for 100/130 mm cryo box	10x10 (vial < 12 mm O.D.)	65 mm	case / 36 pcs..	ML22146

### Cardboard Cryo Boxes 136x136 mm with Dividers

Description	Dimensions	Colour	Sales Unit	Cat. No.
maxxLine® cardboard cryo box B9 incl. 9x9 cells	136x136x50 mm	white	case / 36 pcs..	ML22231
maxxLine® cardboard cryo box B10 incl. 10x10 cells	136x136x50 mm	white	case / 36 pcs..	ML22233

