

ahn myTube®

Screw Cap Tubes



Fully Automated Manufacturing

Manufactured in a controlled Class 100K cleanroom environment with no secondary assembly operation involved reducing the chance of contamination due to manual handling



Unique Infused 'O' Ring Design

Seamlessly moulded during production in Class 100K clean room for reduced contamination risk. Provides superior seal-ability and leakproof stability compared to traditional O-ring designs.

O-ring material: High-quality EPDM (Ethylene Propylene Diene Monomer)



Ergonomic Design

Features an ergonomically knurled surface for a secure grip, ensuring smooth handling even with gloved hands



Built to Endure Rigorous Conditions

Withstands extreme centrifugation speeds up to 20,000g. Durable for -80°C freezing and subsequent thawing at room temperature without compromising integrity



Diverse Volume Ranges & Variants

Available in three convenient volumes
0.5 mL, 1.5 mL, 2 mL

Clear and amber variants, Low retention, skirted or conical formats to suit diverse lab applications



IATA Standards

Tubes are rigorously leak-tested to meet the standards



Sterility Assurance

Screw Cap Tubes are gamma sterilized, ensuring (SAL 10⁻⁶)



Certified

RNase, DNase, Endotoxin,
PCR Inhibitor-free



Cat Code	Description	Bag Qty	Pack size Qty
3-900-50-0	ahn myTube Screw Cap Tube self-standing, 0.5 ml, transparent, Tube without Cap	100	1000 (10x100)
3-901-50-0	ahn myTube Screw Cap Tube self-standing, 0.5 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-902-50-0	ahn myTube Screw Cap Tube self-standing, Low Retention, 0.5 ml, transparent, Tube without Cap	100	1000 (10x100)
3-903-50-0	ahn myTube Screw Cap Tube conical, 1.5 ml, transparent, Tube without Cap	100	1000 (10x100)
3-904-50-0	ahn myTube Screw Cap Tube self-standing, 1.5 ml, transparent, Tube without Cap	100	1000 (10x100)
3-905-50-0	ahn myTube Screw Cap Tube conical, 1.5 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-906-50-0	ahn myTube Screw Cap Tube self-standing, 1.5 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-907-50-0	ahn myTube Screw Cap Tube conical, Low Retention, 1.5 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-908-50-0	ahn myTube Screw Cap Tube self-standing, Low Retention, 1.5 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-909-50-0	ahn myTube Screw Cap Tube conical, 2 ml, transparent, Tube without Cap	100	1000 (10x100)
3-910-50-0	ahn myTube Screw Cap Tube self-standing, 2 ml, transparent, Tube without Cap	100	1000 (10x100)
3-911-50-0	ahn myTube Screw Cap Tube conical, 2 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-912-50-0	ahn myTube Screw Cap Tube self-standing, 2 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-913-50-0	ahn myTube Screw Cap Tube conical, Low Retention, 2 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-914-50-0	ahn myTube Screw Cap Tube self-standing, Low Retention, 2 ml, transparent, sterile, Tube without Cap	100	1000 (10x100)
3-904-50-9	ahn myTube Screw Cap Tube self-standing, 1.5 ml, Amber, Tube without Cap	100	1000 (10x100)
3-906-50-9	ahn myTube Screw Cap Tube self-standing, 1.5 ml, sterile, Amber, Tube without Cap	100	1000 (10x100)
3-908-50-9	ahn myTube Screw Cap Tube self-standing, Low Retention, 1.5 ml, Amber, sterile, Tube without Cap	100	1000 (10x100)

Cat Code	Decription	Bag Qty	Pack size Qty
3-915-50-1	ahn myTube Screw Caps for Screw Cap Tubes, yellow, sterile, Only Caps in a box	100	1000 (10x100)
3-915-50-2	ahn myTube Screw Caps for Screw Cap Tubes, blue, sterile, Only Caps in a box	100	1000 (10x100)
3-915-50-3	ahn myTube Screw Caps for Screw Cap Tubes, red, sterile, Only Caps in a box	100	1000 (10x100)
3-915-50-4	ahn myTube Screw Caps for Screw Cap Tubes, green, sterile, Only Caps in a box	100	1000 (10x100)
3-915-50-9	ahn myTube Screw Caps for Screw Cap Tubes, Amber, sterile, Only Caps in a box	100	1000 (10x100)