

ahn myDPette pro



Bottle Top Dispenser

Premium precision instrument Made in Germany - maximum accuracy and precision, with error ranges better than DIN EN ISO 8655-5 requirements

DE-M marking according to MessEV

Universal dispenser - resistant to most chemicals

Sophisticated system - for air purging and draining without reagent loss

Direct displacement - ideal for crystalline solutions

Built-in recalibration - including a calibration tool

Quick volume adjustment knob - for precise, fast volume setting available for 2.5 mL to 50.0 mL dispenser

2 scales - for measurements in milliliters and ounces

Piston - PTFE protected Viton O-ring displaces liquid during dispensing

Clear view - glass cylinder with protective plastic coating

Screw cap - fully disassemblable for cleaning piston/cylinder

Discharge tube - rotatable with protective sleeve

6 sizes available - from 0.25 mL to 100 mL

Autoclavable at 121°C**



Quick volume adjust knob



Do-it-yourself recalibration

Ordering information

Delivery Package

1 pc. Bottle top dispenser
4 pcs. Thread adapters
2.5/5.0/10.0 mL: A25, A28, A45, S40
25.0/50.0/100.0 mL: A25, A28, A32/45, S40
1 pc. Telescopic suction tube
1 pc. Discharge tube
1 pc. Calibration tool
1 pc. User manual
1 pc. Certificate of performance

Volume	Division	Systematic errors \pm mL	Random errors \leq mL	Thread	Sales Unit	Cat. No.
0.25-2.5 mL	0.05	0.012	0.002	32	carton / 1 pc.	8-710-00-0
0.5-5.0 mL	0.10	0.03	0.005	32	carton / 1 pc.	8-711-00-0
1.0-10.0 mL	0.20	0.06	0.01	32	carton / 1 pc.	8-712-00-0
2.5-25.0 mL	0.50	0.15	0.025	45	carton / 1 pc.	8-713-00-0
5.0-50.0 mL	1.00	0.3	0.05	45	carton / 1 pc.	8-714-00-0
10.0-100.0 mL*	2.00	0.5	0.1	45	carton / 1 pc.	8-715-00-0

* 100 mL version with adjustment screw for volume setting

** read user manual for detailed instructions