

ahn myLab[®]

Single Plate Magnetic Stirrer 1500 RPM SM-01

Usage - in laboratories, pharmaceutical industries, chemical industries

Maximum speed - from 15 to 1500 RPM with an interval of 50 RPM / software-driven gradual acceleration

Time setting - ranging from 1 - 99 minutes / continuous operation feature

Improved mixing - pulse mode where the mixing direction changes every 30 seconds

Maintenance-free design - water spill, corrosion and chemically resistant / includes no moving parts / high-quality stirring bar

User-friendly interface and design - digital display showing the real-time speed and the remaining time / last run memory function / extremely thin (only 16.0 mm thick) / flat surface enables easy wiping and cleaning



Ordering information

Description	Cat. No.	Specifications
myLab [®] SM-01 Single Plate Magnetic Stirrer 1500 RPM	7-300-00-0	Max. stirring quantity 800 mL
		Recommended stir bar length 25 mm
		Protection class IP65
		Dimensions (WxDxH) 156x110x13.2 mm
		Weight 400 g
Optional accessories	Cat. No.	
Magnetic stirring bar - PTFE coated, 25x8 mm	CRS-MAG	
Delivery Package		
1 pc. Magnetic Stirrer 1500 RPM		
1 pc. Stir bar (PTFE coated, 25x8 mm)		
1 pc. Instruction manual		
1 pc. Warranty card		

