


ahn myLab[®]

Tube Roller 30 RPM TR-01



Usage - the mixing of blood samples and viscous substances

Capacity - 9x5 mL tubes

Fixed speed - 30 RPM with an accuracy of ± 2 RPM

Optimal mixing - gentle rocking and rolling movements at a tilt angle of ± 3 degree

Easy maintenance and adjustment - tool free removal and installation of rollers for effortless cleaning and replacing / easy adjustment of spacing to use with different tube diameter

Robust and space-saving design - robust construction for stable operation / small footprint requiring minimum of benchtop space



Tool-free removal and installation of rollers

Ordering information

Description	Cat. No.
myLab [®] TR-01 Tube Roller 30 RPM	7-200-00-0

Specifications	
Run time	Continuous / infinite
Noise level	Noiseless operation
Roller size	30x330 mm
Max. loading capacity	1.5 kg
Protection class	IP 21
Dimensions (WxDxH)	475x201x80 mm
Weight	2.0 kg

Delivery Package
1 pc. Tube Roller 30 RPM
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card