

ahn myPlate®



Deep Well Plates

1.0 mL | 2.0 mL | 2.2 mL

Contaminant-free - manufactured in a class 100K cleanroom facility

Certified purity - free of DNase, RNase and endotoxin inhibitors

Fully autoclavable - thermostable to 121°C, AHN myPlate® Deep Well Plates are autoclave-ready for easy sterilisation and decontamination after use. Safe for use in high-temperature protocols up to 121°C

Freezer-ready - excellent low temperature stability, AHN myPlate® Deep Well Plates can be used for sample and reagent storage in refrigerators and freezers. Safe for use on low temperatures up to a limit of -80°C


Alphanumeric sample index - feature an intuitive, embossed sample index with alphanumeric code for easy well identification

ANSI compliant - including ANSI/SLAS standards for microplate dimensions ensuring universal compatibility with standard handling platforms

Easy stacking - compact design saves space on workbench or in storage

Raised-rim design - minimises cross-contamination across sample wells, and feature a subtle raised-rim design around the opening of each well

Compatibility list - see page 97



Ordering information

Volume	Bottom	Top	Colour	Sales Unit	Cat. No.
myPlate® Deep Well Plate 1.0 mL, Round bottom		round	clear	case / 4x24	3-126-01-0
myPlate® Deep Well Plate 1.0 mL, Round bottom NEW PRODUCT		square	clear	case / 4x24	3-126-02-0
myPlate® Deep Well Plate 2.0 mL, Round bottom		round	clear	case / 4x12	3-226-01-0
myPlate® Deep Well Plate 2.0 mL, V-bottom		square	clear	case / 4x12	3-226-02-0
myPlate® Deep Well Plate 2.2 mL, Zybio®, U-bottom		square	clear	case / 10x5	3-226-03-0
myPlate® Deep Well Plate 2.2 mL, Kingfisher®, V-bottom		square	clear	case / 10x5	3-226-04-0
myPlate® 96 Kingfisher® type tip comb, V-bottom		round	clear	case / 100x2	3-230-00-0
myPlate® 8-well Zybio® type tip comb, U-bottom		round	clear	case / 100x2	3-230-01-0

Kingfisher is a trademark of Thermo Fisher Scientific Inc. Zybio is a trademark of Zybio Inc.

ahn myPlate

Silicone mat

Ordering information

Volume	Sales Unit	Cat. No.
myPlate® Silicone mat All fit for square top cross cut	case / 10x50 pcs.	3-231-00-0
myPlate® Silicone mat Round 1.0 mL pierceable	case / 10x50 pcs.	3-231-01-0
myPlate® Silicone mat Round 2.0 mL pierceable	case / 10x50 pcs.	3-231-02-0

