



**Superior quality, contaminant-free petri dishes** - manufactured in a cleanroom facility, resulting in sterile dishes

**Certified purity** - free of endotoxin and pyrogen

**Space-saving** - with a standard sized diameter of 90 mm, AHN Petri Dishes have the perfect balance between plating area and space-saving dimensions. These dishes benefit from being both lightweight and fully stackable - raised rim allows easy stacking by hand without slipping or scratching

**Excellent airflow** - supplied in vented and non-vented formats to accommodate different microbial culture needs. Vented plates have a lip that allows for airflow and moisture control whilst non-vented dishes limit airflow and maximise moisture retention.

**High-quality material** - made of high-quality polystyrene that meets the specifications of US FDA regulation 21 CFR 177.1640

**Product traceability** - each batch is supplied with a traceable lot number as well as information on sterility and expiration

**Compatibility** - fully compatible with medium-throughput carousel-based media preparation systems / Compatibility Chart is available on request

**Packaging Options** - available in 24x20 packs with grid and 25x20 packs



## Ordering information

Product	Diameter	Height	Vented	Sterile	Colour	Pack Type	Sales Unit	Cat. No.
myPlate® Petri Dish	91.5 mm	14.35 mm	no	■	clear	bag	case / 24x20	3-520-48-0
myPlate® Petri Dish	91.5 mm	14.35 mm	no	■	clear	bag	case / 24x20	3-521-48-0
myPlate® Petri Dish	91.5 mm	14.55 mm	yes, 3 cams	aseptic	clear	bag	case / 24x20	3-522-48-0
myPlate® Petri Dish	91.5 mm	14.55 mm	yes, 3 cams	■	clear	bag	case / 24x20	3-523-48-0
myPlate® Petri Dish	91.5 mm	14.35 mm	no	aseptic	clear	bag	case / 25x20	3-520-50-0
myPlate® Petri Dish	91.5 mm	14.35 mm	no	■	clear	bag	case / 25x20	3-521-50-0
myPlate® Petri Dish	91.5 mm	14.55 mm	yes, 3 cams	aseptic	clear	bag	case / 25x20	3-522-50-0
myPlate® Petri Dish	91.5 mm	14.55 mm	w/ yes, 3 cams	■	clear	bag	case / 25x20	3-523-50-0