


ahn myTip®

Effortless EL³ Tips

AHN myTip Effortless EL³ Tips are designed for optimal performance with Rainin® LTS® style pipettes. The innovative design of the tips is characterised by a cylindrical collar with a positive stop feature, which enables the tips to be attached and removed quickly and effortlessly. At the same time, it ensures a perfect seal between the pipette and the tip.

AHN myTip Effortless EL³ Tips provide the accuracy, precision and ergonomics you expect when using your Rainin® LTS® style pipette.



Designed for optimal performance with Rainin® LTS® style pipettes



Certified free of RNase, DNase, PCR Inhibitors, Endotoxins and Pyrogens



Cylindrical shape and positive stop ensure a perfect fit

US FDA approved medical grade virgin polypropylene plastic without any trace of metals, highly transparent and clear



Colour-coded trays for easy identification of tip and pipette volume

Class 100000 clean room facility with automated manufacturing



Graduation marks used for visual control to minimise the risk of errors during pipetting



Available in rack, reload and Eco-rack

ahn

Uthleber Weg 14
99734 Nordhausen
www.ahn-bio.de

P: +49(0)3631/65242-0
F: +49(0)3631/65242-90
E: info@ahn-bio.com

AHN myTip Effortless EL³ - Tips



20 µL

Cat No.	Description	Sterile	Sales Unit
1-711-C1-R	AHN myTip Effortless EL ³ Tips 20 µL, reload		Case / 10 x 10 x 96 tips
1-712-C1-R	AHN myTip Effortless EL ³ Tips 20 µL, reload	■	Case / 10 x 10 x 96 tips
1-716-C10-0	AHN myTip Effortless EL ³ Tips 20 µL, racked		Case / 10 x 10 x 96 tips
1-718-C10-0	AHN myTip Effortless EL ³ Tips 20 µL, racked	■	Case / 10 x 10 x 96 tips
1-716-E10-0	AHN myTip Effortless EL ³ Tips 20 µL, Eco-rack		Case / 10 x 10 x 96 tips
1-718-E10-0	AHN myTip Effortless EL ³ Tips 20 µL, Eco-rack	■	Case / 10 x 10 x 96 tips



200 µL

Cat No.	Description	Sterile	Sales Unit
1-721-C10-0	AHN myTip Effortless EL ³ Tips 200 µL, reload		Case / 10 x 10 x 96 tips
1-722-C10-0	AHN myTip Effortless EL ³ Tips 200 µL, reload	■	Case / 10 x 10 x 96 tips
1-721-C1-R	AHN myTip Effortless EL ³ Tips 200 µL, racked		Case / 10 x 10 x 96 tips
1-722-C1-R	AHN myTip Effortless EL ³ Tips 200 µL, racked	■	Case / 10 x 10 x 96 tips
1-721-E10-0	AHN myTip Effortless EL ³ Tips 200 µL, Eco-rack		Case / 10 x 10 x 96 tips
1-722-E10-0	AHN myTip Effortless EL ³ Tips 200 µL, Eco-rack	■	Case / 10 x 10 x 96 tips



1000 µL

Cat No.	Description	Sterile	Sales Unit
1-802-C10-0	AHN myTip Effortless EL ³ Tips 1000 µL, reload		Case / 10 x 6 x 96 tips
1-803-C10-0	AHN myTip Effortless EL ³ Tips 1000 µL, reload	■	Case / 10 x 6 x 96 tips
1-802-C4-R	AHN myTip Effortless EL ³ Tips 1000 µL, racked		Case / 10 x 6 x 96 tips
1-803-C4-R	AHN myTip Effortless EL ³ Tips 1000 µL, racked	■	Case / 10 x 6 x 96 tips
1-802-E10-0	AHN myTip Effortless EL ³ Tips 1000 µL, Eco-rack		Case / 10 x 6 x 96 tips
1-803-E10-0	AHN myTip Effortless EL ³ Tips 1000 µL, Eco-rack	■	Case / 10 x 6 x 96 tips

AHN myTip Effortless EL³ - Filter Tips



20 µL

Cat No.	Description	Sterile	Sales Unit
2-718-C10-0	AHN myTip Effortless EL ³ Filter Tips, 20 µL, racked	■	Case / 10 x 10 x 96 tips
2-718-E10-0	AHN myTip Effortless EL ³ Filter Tips, 20 µL, Eco-rack	■	Case / 10 x 10 x 96 tips



200 µL

Cat No.	Description	Sterile	Sales Unit
2-728-C10-0	AHN myTip Effortless EL ³ Filter Tips, 200 µL, racked	■	Case / 10 x 10 x 96 tips
2-728-E10-0	AHN myTip Effortless EL ³ Filter Tips, 200 µL, Eco-rack	■	Case / 10 x 10 x 96 tips



1000 µL

Cat No.	Description	Sterile	Sales Unit
2-803-C10-0	AHN myTip Effortless EL ³ Filter Tips, 1000 µL, racked	■	Case / 10 x 6 x 96 tips
2-803-E10-0	AHN myTip Effortless EL ³ Filter Tips, 1000 µL, Eco-rack	■	Case / 10 x 6 x 96 tips