


ahn *myVial*[®]

Diagnostic Vial



Fully Automated Manufacturing

Manufactured from 100K class cleanroom
Constructed from USP Class VI US FDA
Approved PETG material



Exceptional Clarity

Ultra-clear PETG ensures effortless visual
inspection of contents, making it ideal for
diagnostic and clinical applications.



Diverse Volume Ranges

Available in three volumes 5ml, 10ml, 20ml



Leakproof Security

Each vial is fitted with a precision-engineered
HDPE screw cap and an integrated polyethylene
cap liner that ensures a secure seal, minimizing
leakage during transport.



Non-Cytotoxic & Non-Pyrogenic

Tested and certified to be safe for biological
samples, supporting sensitive diagnostic and
pharmaceutical workflows.



Packaging

Available in a variety of standard sizes with
tailored packaging options to match your
workflow needs



Sterility Assurance

Gamma sterilized (SAL 10⁻⁶)
in ISO Class 7 cleanroom
conditions for assured
contamination-free
performance.



Certified

RNase, DNase, Endotoxin,
PCR Inhibitor-free



Ordering information

Cat Code	Decription	Quantity/Case
AHN-PETG-10-S	ahn® myVial Diagnostic Vial PETG 10 ml, sterile, racked	5 x 20 pcs
AHN-PETG-10-S-B	ahn® myVial Diagnostic Vial PETG 10 ml, sterile, racked, bulk	5 x 100 pcs
AHN-PETG-20-S	ahn® myVial Diagnostic Vial PETG 20 ml, sterile, racked	5 x 20 pcs
AHN-PETG-20-S-B	ahn® myVial Diagnostic Vial PETG 20 ml, sterile, racked, bulk	5 x 100 pcs
AHN-PETG-5-S	ahn® myVial Diagnostic Vial PETG 5 ml, sterile, racked	5 x 20 pcs
AHN-PETG-5-S-B	ahn® myVial Diagnostic Vial PETG 5 ml, sterile, racked, bulk	5 x 100 pcs