

ahn pipet4u[®]

Oasis

Pipette Controller

UV
RESISTANT

2
YEARS
WARRANTY

autoclavable
121°C
15 min

Comfortable and effortless pipetting - lightweight pipette design / smooth push buttons for precise pipetting control / ergonomic design for easy operation

Cordless working - holds charge up to 8 hours / quiet high-powered motor

Contamination control - replaceable membrane filter to minimise risk of contamination

Speed setting - adjustment of pipetting speed on sample viscosity with three selectable speeds

Sterile handling - autoclavable pipette cone (121°C / 15 min)



Autoclavable tip cone



Ergonomic push buttons

Ordering information

Description	Sales Unit	Cat. No.
pipet4u [®] Oasis Pipette Controller 0.1-100 mL	carton / 1 pc.	8-300-00-9

Description - Spare Parts	Sales Unit	Cat. No.
0.20 µm PTFE membrane filter	bag / 5 pcs.	8-300-03-0
0.45 µm PTFE membrane filter	bag / 5 pcs.	8-300-04-0
Silicone adapter	bag / 1 pc.	8-300-01-0
Pipette cone	bag / 1 pc.	8-300-02-8

Delivery Package
1 pc. Pipette Controller
2 pcs. PTFE Membrane filters (0.2 µm / 0.45 µm)
1 pc. Charger (EU)
2 pcs. NiMH batteries
1 pc. Wall mount bracket
1 pc. Desktop stand
1 pc. Instruction manual
1 pc. Warranty card



+49(0)3631/65242-0



contact@ahn-bio.com



www.ahn-bio.com

