

ahn *reagent* Reservoirs

Construction - formulated from of a unique blend of polystyrene

Extra large, sturdy base with broad edges - enables it to stand firm on your work bench

Graduation marks - helps identify any remaining liquid and facilitates easy, accurate dispensing

Angled interior - helps the collection of the liquid and decreases reagent waste

Spillproof - no need to worry about sample spillage while you are processing liquids

Easy pouring - to store any remaining liquid you can empty the reservoir using the corner spouts

Compact storage - easy to stack and the lip at the bottom of the reservoir allows them to be easily separated from each other






Ideal for using with all AHN Single and Multi-channel Pipettes

- our reservoirs come in three different sizes, to suit various pipettes:

- ① 10 mL, V-shaped bottom (ideal with 8- and 16-channel pipettes)
- ② 25 mL, flat bottom (ideal with multi-channel pipettes)
- ③ 50 mL, V-shaped bottom (ideal with 12-channel pipettes)



Ordering information

| Description | Colour | Pack Type | Sales Unit | Cat. No. |
|---|--------|-----------|--------------|------------|
| Reagent Reservoirs 10 mL  | white | bag | case / 5x100 | 3-701-11-9 |
| Reagent Reservoirs 25 mL  | white | bag | case / 5x100 | 3-702-11-9 |
| Reagent Reservoirs 50 mL  | white | bag | case / 5x100 | 3-703-11-9 |



+49(0)3631/65242-0



contact@ahn-bio.com



www.ahn-bio.com

