

ahn[®]myPette[®] EL³

Manual Pipettes

EL³ technology uses cylindrical tip and shaft mechanism - eliminates high force to deliver much more comfortable pipetting experience

Effortless pipetting - optimal weight distribution and lowest plunger force requirement through innovative spring and seal design

Effortless tip ejection - to ensure continuous uninterrupted pipetting without any strain tip ejection

Magnet-assisted piston - consistently precise results

Effortless tip loading - for a reliable and consistent seal between the pipette and the tip without applying excessive force

Colour coding - for easy identification 

Volume lock - prevents any accidental changes in the volume and guarantees constant pipetting


Large 4-digit display - Precise setting of the volume through an easy-to-read display

Highly durable material - to protect against chemical corrosion and physical stresses

Calibrated - as per ISO 8655-6

Fully autoclavable and UV resistant

**NEW
PRODUCT**



Delivery Package

1 pc. Manual pipette	1 pc. Shelf-mounting stand
1 pc. Conformity/calibration certificate	1 pc. Calibration tool
1 pc. Warranty card	1 pc. Silicone grease
1 pc. Instruction manual	1 pc. Sample pack of tips

¹ Systematic error² Random error

Ordering information

Single Channel Pipette, Variable Volume

Volume	Colour Code	Channels	Test volume	Increments	Inaccuracy ¹	Imprecision ²	Cat. No.
0.1-2.5 µL		1	0.25	0.002	12.00	6.00	8-100-32-9
			1.25		2.50	1.20	
			2.5		1.50	0.69	
0.5-10 µL		1	1	0.02	2.50	1.20	8-101-32-9
			5		1.50	0.60	
			10		1.00	0.40	
2-20 µL		1	2	0.02	7.00	2.00	8-102-32-9
			10		1.50	0.50	
			20		1.00	0.30	
5-50 µL		1	5	0.1	3.50	1.50	8-103-32-9
			25		1.20	0.40	
			50		0.80	0.20	
10-100 µL		1	10	0.1	3.50	1.00	8-104-32-9
			50		0.80	0.21	
			100		0.70	0.15	
20-200 µL		1	20	0.2	2.50	1.00	8-105-32-9
			100		0.70	0.23	
			200		0.70	0.15	

ahn[®]



Uthleber Weg 14 | 99734 Nordhausen | Germany

Single Channel Pipette, Variable Volume

Volume	Colour Code	Channels	Test volume	Increments	Inaccuracy ¹	Imprecision ²	Cat. No.
100-1000 µL		1	100	1	2.50	0.60	8-106-32-9
			500		0.80	0.20	
			1000		0.80	0.15	

Multi Channel Pipette, Variable Volume

Volume	Colour Code	Channels	Test volume	Increments	Inaccuracy ¹	Imprecision ²	Cat. No.
0.5-10 µL		8	1	0.02	2.50	1.20	8-201-32-9
			5		1.50	0.60	
			10		1.00	0.40	
2-20 µL		8	20	0.02	2.50	1.00	8-205-32-9
			100		0.70	0.23	
			200		0.70	0.15	
5-50 µL		8	2	0.1	7.00	2.00	8-202-32-9
			10		1.50	0.50	
			20		1.00	0.30	
10-100 µL		8	10	0.1	3.50	1.00	8-204-32-9
			50		0.80	0.21	
			100		0.70	0.15	
20-200 µL		8	30	0.2	2.50	1.00	8-206-32-9
			150		0.80	0.25	
			300		0.80	0.15	
30-300 µL		8	5	0.2	3.50	1.50	8-203-32-9
			25		1.20	0.40	
			50		0.80	0.20	
0.5-10 µL		12	1	0.02	2.50	1.20	8-401-32-9
			5		1.50	0.60	
			10		1.00	0.40	
2-20 µL		12	20	0.02	2.50	1.00	8-405-32-9
			100		0.70	0.23	
			200		0.70	0.15	
5-50 µL		12	2	0.1	7.00	2.00	8-402-32-9
			10		1.50	0.50	
			20		1.00	0.30	
10-100 µL		12	10	0.1	3.50	1.00	8-404-32-9
			50		0.80	0.21	
			100		0.70	0.15	
20-200 µL		12	30	0.2	2.50	1.00	8-406-32-9
			150		0.80	0.25	
			300		0.80	0.15	
30-300 µL		12	5	0.2	3.50	1.50	8-403-32-9
			25		1.20	0.40	
			50		0.80	0.20	

¹ Systematic error² Random error



+49(0)3631/65242-0



contact@ahn-bio.com



www.ahn-bio.com

