

ahn® myReservoir®

Foldable Reagent Reservoirs

30 mL/40 mL

V-shaped bottom - prevents bubble formation and ensures retrieval of the last drop of liquid for efficient reagent use

Sturdy construction - offers a consistently firm and stable reservoir for pipetting

Pre-folded design - allows easy assembly by simply pushing the sides together

Materials - Made from durable polypropylene, ensuring reliability and reusability


Cost-saving - Shipped and stored as flat sheets, saving space and reducing shipping costs

Available Volumes - as 30 mL reservoirs for 8/16 channel pipettes and 40 mL reservoirs for 12 channel pipettes


Eco-friendly - and space-saving, optimising storage and reducing environmental impact

Space- Saving Increases - Increases workspace efficiency by saving benchtop space and lowering costs


Ready in three easy steps...



Remove a flat reservoir from pack



Push two outer sides inwards towards the centre



Push end triangles towards centre to complete reservoir assembly



49

Ordering information

Description	Sterile	Volume	For use with	Sales Unit	Cat. No.
ahn® myReservoir Foldable Reagent Reservoir	■	30 mL	8/16 channel pipettes	case / 10x5	3-711-26-0
ahn® myReservoir Foldable Reagent Reservoir	■	40 mL	12 channel pipettes	case / 10x5	3-713-26-0



+49(0)3631/65242-0



contact@ahn-bio.com



www.ahn-bio.com

