

Designed for optimal performance - with Rainin[®] LTS[®] style pipettes

Cylindrical tip and shaft mechanism - innovative tip design of the tips ensures a perfect fit, and EL³ technology eliminates high force to deliver much more comfortable pipetting experience

Positive stop feature enables the tips to be attached and removed quickly and effortlessly, while ensuring a perfect seal between the pipette and the tip



Colour-coded trays - for easy identification of tip and pipette volume

Graduation marks used for visual control - to facilitate identification and minimise the risk of errors during pipetting

Materials - Certified free of RNase, DNase, PCR Inhibitors, Endotoxins and Pyrogens. US FDA approved medical grade virgin polypropylene plastic without any trace of metals, highly transparent and clear

Manufactured - in Class ISO 8 cleanroom facility with automated manufacturing

Versions - available in rack, reload, eco-rack








Ordering information

Pipette Tips













Volume	Length	Colour Tip	Colour Code	Sterile	Pack Type	Sales Unit	Cat No.
 myTip [®] EL ³ tips, 20 µL	49.4 mm	clear	■	■	reload	case / 10x10x96	1-711-C1-R
				■	reload	case / 10x10x96	1-712-C1-R
				■	racked	case / 10x10x96	1-716-C10-0
				■	eco-rack	case / 10x10x96	1-716-E10-0
				■	racked	case / 10x10x96	1-718-C10-0
				■	eco-rack	case / 10x10x96	1-718-E10-0
 myTip [®] EL ³ tips, 200 µL	50.0 mm	clear	■	■	racked	case / 10x10x96	1-721-C10-0
				■	reload	case / 10x10x96	1-721-C1-R
				■	eco-rack	case / 10x10x96	1-721-E10-0
				■	racked	case / 10x10x96	1-722-C10-0
				■	reload	case / 10x10x96	1-722-C1-R
				■	eco-rack	case / 10x10x96	1-722-E10-0

LTS[®] is a registered trademark of Rainin Instrument. Rainin[®] is a registered trademark of Mettler Toledo Rainin LLC
















Volume	Length	Colour Tip	Colour Code	Sterile	Pack Type	Sales Unit	Cat No.
 myTip® EL ³ tips, 1000 µL	85.8 mm	clear			racked	case / 10x6x96	1-802-C10-0
					reload	case / 10x6x96	1-802-C4-R
					eco-rack	case / 10x6x96	1-802-E10-0
					racked	case / 10x6x96	1-803-C10-0
					reload	case / 10x6x96	1-803-C4-R
					eco-rack	case / 10x6x96	1-803-E10-0













Filter Tips

Volume	Length	Colour Tip	Colour Code	Sterile	Pack Type	Sales Unit	Cat No.
 myTip® EL ³ filter tips, 20 µL	49.4 mm	clear			racked	case / 10x10x96	2-718-C10-0
					eco-rack	case / 10x10x96	2-718-E10-0
 myTip® EL ³ filter tips, 200 µL	50.0 mm	clear			racked	case / 10x10x96	2-728-C10-0
					eco-rack	case / 10x10x96	2-728-E10-0
 myTip® EL ³ filter tips, 1000 µL	85.8 mm	clear			racked	case / 10x6x96	2-803-C10-0
					eco-rack	case / 10x6x96	2-803-E10-0

Low Retention Tips

Description	Length	Colour Tip	Colour Code	Sterile	Pack type	Sales Unit	Cat No.
 myTip® EL ³ tips, 20 µL, low retention	49.4 mm	clear			reload	case / 10x10x96	4-711-C1-R
					racked	case / 10x10x96	4-718-C10-0
 myTip® EL ³ tips, 200 µL, low retention	50.0 mm	clear			reload	case / 10x10x96	4-721-C1-R
					racked	case / 10x10x96	4-728-C10-0
					eco-rack	case / 10x10x96	4-728-E10-0
 myTip® EL ³ tips, 1000 µL, low retention	85.8 mm	clear			reload	case / 10x10x96	4-802-C4-R
					racked	case / 10x10x96	4-803-C10-0
					eco-rack	case / 10x10x96	4-803-E10-0

Low Retention Filter Tips

Description	Length	Colour Tip	Colour Code	Sterile	Pack type	Sales Unit	Cat No.
 myTip® EL ³ filter tips, 20 µL, low retention	49.4 mm	clear			racked	case / 10x10x96	5-718-C10-0
					eco-rack	case / 10x10x96	5-718-E10-0
 myTip® EL ³ Filter Tips, 200 µL, low retention	50.0 mm	clear			racked	case / 10x10x96	5-728-C10-0
					eco-rack	case / 10x10x96	5-728-E10-0
 myTip® EL ³ Filter Tips, 1000 µL, low retention	85.8 mm	clear			racked	case / 10x6x96	5-803-C10-0
					eco-rack	case / 10x6x96	5-803-E10-0

LTS® is a registered trademark of Rainin Instrument. Rainin® is a registered trademark of Mettler Toledo Rainin LLC