


# CAPP *Bravo*

SINGLE | MULTICHANNEL PIPETTES

CAPP Bravo pipettes combine lightweight ergonomic design, and the utmost accuracy and precision in an affordable package

## Features

- Ultra lightweight construction
- Metallic piston results in accurate liquid dispensing, delivering consistent, reliable results across all channels
- Built-in oversized volume indicator simplifies volume monitoring and makes precise adjustment particularly easy
- Durable construction ensures chemical resistance, thermal stability and resilience from accidental bumps or damage
- Long-term calibration retention, even after sterilising in an autoclave
- Metal shaft with double O-rings for all multichannel volumes and 0.1-2  $\mu$ L, 0.5-10  $\mu$ L, 2-20  $\mu$ L single channel volumes
- Coloured caps and identification stickers enable the easy visual identification of pipettes or volume capacities



**2 YEAR**  
Warranty

**Fully**  
Autoclavable

# Ordering Information

## CAPPBravo Single Channel Pipettes - Variable Volume

| Cat. No. | Channels | Description | Increments | Inaccuracy % | Imprecision % | Colour                                                                              |
|----------|----------|-------------|------------|--------------|---------------|-------------------------------------------------------------------------------------|
| B02-1    | 1        | 0.1-2 µL    | 0.01µl     | 12.00 / 2.00 | 6.00 / 1.00   |  |
| B10-1    | 1        | 0.5-10 µL   | 0.02µl     | 2.50 / 1.00  | 1.80 / 0.50   |  |
| B20-1    | 1        | 2-20 µL     | 0.02µl     | 5.00 / 1.00  | 1.50 / 0.30   |  |
| B50-1    | 1        | 5-50 µL     | 0.1µl      | 2.00 / 1.00  | 1.20 / 0.30   |  |
| B100-1   | 1        | 10-100 µL   | 0.1µl      | 2.50 / 0.70  | 1.00 / 0.20   |  |
| B200-1   | 1        | 20-200 µL   | 0.2µl      | 2.50 / 0.60  | 0.70 / 0.20   |  |
| B1000-1  | 1        | 100-1000 µL | 1.0µl      | 2.00 / 0.80  | 0.60 / 0.20   |  |
| B5000-1  | 1        | 1-5 mL      | 10.0µl     | 1.20 / 0.60  | 0.30 / 0.15   |  |
| B10000-1 | 1        | 1-10 mL     | 20.0µl     | 3.00 / 0.60  | 0.60 / 0.15   |  |

## CAPPBravo Single Channel Pipettes - Fixed Volume

| Cat. No.  | Channels | Description | Increments | Inaccuracy % | Imprecision % | Colour                                                                                |
|-----------|----------|-------------|------------|--------------|---------------|---------------------------------------------------------------------------------------|
| B01-1F    | 1        | 1 µL        | -          | 2.50         | 1.80          |    |
| B02-1F    | 1        | 2 µL        | -          | 2.00         | 1.00          |    |
| B05-1F    | 1        | 5 µL        | -          | 1.20         | 0.80          |    |
| B10-1F    | 1        | 10 µL       | -          | 1.00         | 0.50          |   |
| B20-1F    | 1        | 20 µL       | -          | 1.00         | 0.30          |  |
| B50-1F    | 1        | 50 µL       | -          | 1.00         | 0.30          |  |
| B100-1F   | 1        | 100 µL      | -          | 0.70         | 0.20          |  |
| B200-1F   | 1        | 200 µL      | -          | 0.60         | 0.20          |  |
| B500-1F   | 1        | 500 µL      | -          | 0.80         | 0.15          |  |
| B1000-1F  | 1        | 1000 µL     | -          | 0.80         | 0.20          |  |
| B5000-1F  | 1        | 5 mL        | -          | 0.60         | 0.15          |  |
| B10000-1F | 1        | 10 mL       | -          | 0.60         | 0.15          |  |

## CAPPBravo Multichannel Pipettes

| Cat. No. | Channels | Description | Increments | Inaccuracy % | Imprecision % | Colour                                                                                |
|----------|----------|-------------|------------|--------------|---------------|---------------------------------------------------------------------------------------|
| B10-8    | 8        | 0.5-10 µL   | 0.02 µl    | 5.00 / 2.00  | 3.60 / 1.00   |  |
| B100-8   | 8        | 10-100 µL   | 0.1 µl     | 2.50 / 1.30  | 1.20 / 0.50   |  |
| B300-8   | 8        | 30-300 µL   | 0.2 µL     | 1.66 / 1.40  | 0.60 / 0.30   |  |
| B10-12   | 12       | 0.5-10 µL   | 0.02 µl    | 5.00 / 2.00  | 3.60 / 1.00   |  |
| B100-12  | 12       | 10-100 µL   | 0.1 µl     | 2.50 / 1.30  | 1.20 / 0.50   |  |
| B300-12  | 12       | 30-300 µL   | 0.2 µL     | 1.66 / 1.40  | 0.60 / 0.30   |  |

## CAPPBravo Starter Kits

| Cat. No.  | Description                                                                          |
|-----------|--------------------------------------------------------------------------------------|
| SK-03-B01 | Starter Kit including B10-1 / B100-1 / B1000-1 pipettes, C-05 stand, 3 racks of tips |
| SK-03-B02 | Starter Kit including B20-1/B200-1/B1000-1 pipettes, C-05 stand, 3 racks of tips     |
| SK-03-BB  | Three pipettes free of choice, C-05 stand, 3 racks of tips                           |