



EXPell

CRYO TUBES

CAPP offers EXPell Cryo Tubes in 5 sizes: 0.5 mL, 1.0 mL, 2.0 mL, 4.0 mL and 5.0 mL. Focusing on the best user experience and convenience, the quick seal closure with external threads enables easy tightening or removal with a half-turn.

These are one of very few cryo tubes to have external threads and star feet, allowing storage in 10x10 cryo boxes and saving up to 23% space compared to other standard options on the market. Made from specially formulated polypropylene, both the caps and tubes have matching expansion coefficients, ensuring leakproof performance at various temperatures.

EXPell Cryo Tubes are ideal for storage in the vapour phase of liquid nitrogen.

Coloured cap inserts. Available in red, blue, green, white and yellow. Made of PP for easy recognition

Lip seal cap design ensures that the tubes are **leakproof**

High-quality graduation marks and writing area for **easy sample identification**

Super straight design fits into both 9x9 and 10x10 cryo boxes. No surface transition between the tube and the cap **for minimal risk of contamination during storage**



External thread super straight design gives **up to 23% extra storage capacity** when used with the 10x10 cryo boxes

Certified Purity – free of RNase, DNase, human DNA, PCR-Inhibitors, Pyrogens and ATP

Quick seal design for single-handed operation. Tighten or open it with only half a turn!

Star-foot bottom for better stability and operation when stored on selected cryo boxes

Ordering Information

EXPell Cryo Tubes

Cat. No.	Description
5040500C	EXPell Cryo Tube 0,5 mL, sterile, w/o print, bulk, case with 10x100 pcs (tubes) and 1x1000 pcs (caps)
5041000C	EXPell Cryo Tube 1.0 mL, 10x100 pcs.
5042000C	EXPell Cryo Tube 2.0 mL, 10x100 pcs.
5044000C	EXPell Cryo Tube 4.0 mL, 5x100 pcs.
5045000C	EXPell Cryo Tube 5.0 mL, 5x100 pcs.

Coloured Cap Inserts for EXPell Cryo Tubes

Cat. No.	Description
504R000C	Red Cap Insert for EXPell Cryo Tube, 1x500 pcs.
504B000C	Blue Cap Insert for EXPell Cryo Tube, 1x500 pcs.
504G000C	Green Cap Insert for EXPell Cryo Tube, 1x500 pcs.
504W000C	White Cap Insert for EXPell Cryo Tube, 1x500 pcs.
504Y000C	Yellow Cap Insert for EXPell Cryo Tube, 1x500 pcs.