



# EXPell EXPell *Secure*

MICROCENTRIFUGE TUBES

With EXPell and EXPellSecure, CAPP offers two variants of microcentrifuge tubes. Both manufactured from high-quality raw materials that meet standards of purity and have high mechanical stability and excellent chemical resistance. The EXPellSecure microcentrifuge tubes feature an exclusive design that eliminates the risk of accidental opening or lid detachment during vigorous processes. Hinge orientation nibs and safe click lock ensure leakproof performance and reducing sample loss due to evaporation during incubation or storage.

## Features

- Perfect seal prevents evaporation and ensures sample integrity
- Easy opening/closing with one hand
- Quick volume control thanks to precise graduation marks
- Frosted areas on tube body and lid allow easy sample identification
- For colour coding of the samples, the 1.5 mL EXPellSecure microcentrifuge tubes are also available in colour-coded sales units
- High clarity for excellent sample visibility
- Free of DNase, RNase, human DNA, pyrogen and endotoxin
- For high-performance applications EXPell Microcentrifuge Tubes withstands up to 20000 x g and EXPellSecure Microcentrifuge Tube up to 25000 x g



## Ordering Information

Cat. No.	Description
5100500C	EXPell Microcentrifuge tubes 0.5 mL, clear, bulk, 10x1000 pcs.
5100505C	EXPell Microcentrifuge tubes 0.5 mL, clear, bulk, sterile, 20x500 pcs.
5101500C	EXPellSecure Microcentrifuge tubes 1.5 mL, clear, bulk, 16x500 pcs.
5101501C	EXPell Microcentrifuge tubes 1.5 mL, clear, bulk, 20x500 pcs.
5101505C	EXPell Secure Microcentrifuge tubes 1.5 mL, clear, sterile, bulk, 10x500 pcs.
5101590C	Expell Secure microcentrifuge tubes 1.5 mL, assorted colours, bulk, 20x500 pcs. (magenta, blue, yellow, green, orange)
5102000C	EXPellSecure Microcentrifuge tubes 2.0 mL, clear, bulk, 16x500 pcs.
5102001C	EXPell Microcentrifuge tubes 2.0 mL, clear, bulk, 10x500 pcs.
5102005C	EXPell Secure microcentrifuge tubes 2.0 mL, clear, sterile, bulk, 10x500 pcs.