


CAPP Origami

REAGENT RESERVOIRS



CAPPOrigami represents a new way of thinking when it comes to reagent reservoirs.

Features

- V-shaped bottom prevents bubble formation and ensures retrieval of the last drop of liquid for efficient reagent use
- Sturdy construction offers a consistently firm and stable reservoir for pipetting
- Pre-folded design allows easy assembly by simply pushing the sides together
- Made from durable polypropylene, ensuring reliability and reusability
- Shipped and stored as flat sheets, saving space and reducing shipping costs
- Reservoirs are available as 30 mL reservoirs for 8/16 channel pipettes and 40 mL reservoirs for 12 channel pipettes
- Eco-friendly and space-saving, optimising storage and reducing environmental impact
- Increases workspace efficiency by saving benchtop space and lowering costs



Ordering Information

Cat. No.	Description
CA40505	CAPPOrigami 30 mL (8- and 16-channel pipettes), bag w/ 50 pcs.
CA40506	CAPPOrigami 30 mL (8- and 16-channel pipettes), sterile, 10 bags w/ 5 pcs. each
CA40510	CAPPOrigami 40 mL (12-channel pipettes), bag w/ 50 pcs.
CA40511	CAPPOrigami 40 mL (12-channel pipettes), sterile, 10 bags w/ 5 pcs. each