

CAPP offers a comprehensive selection of practical accessories to complement and enhance the functionality of its full line of pipettes.

Volume knobs


- **Pre-calibrated interchangeable**, variable and fixed volume knobs
- Allow users to **use one pipette with many fixed volume knobs** to substantially reduce cost. Now, instead of four fixed volume pipettes, for example, user can buy one CAPP pipette and four inexpensive knobs to achieve the same precise results.

Pipette Stands

- Delivers optimal storage solutions for laboratories
- Options hold between **3 and 6 mechanical pipettes**

Pipette Filters

- An economical solution to eliminate “splash effects” and will substantially reduce risk of sample and shaft contamination



Ordering Information

Fixed Volume Knobs - for CAPP ecopipette, CAPP Aero 96 and CAPP Aero 384

Cat. No.	Volume	Use with pipette volumes	Colour
02-FE-1	1 µL	0.2-2 µL	●
02-FE-2	2 µL	0.2-2 µL	●
10-FE-5	5 µL	0.5-10 µL	○
10-FE-10	10 µL	0.5-10 µL	○
20-FE-10	10 µL	2-20 µL	●
20-FE-20	20 µL	2-20 µL	●
50-FE-20	20 µL	5-50 µL	●
50-FE-25	25 µL	5-50 µL	●
50-FE-50	50 µL	5-50 µL	●
100-FE-30	30 µL	10-100 µL	●
100-FE-50	50 µL	10-100 µL	●
100-FE-100	100 µL	10-100 µL	●

Fixed Volume Knobs - for CAPP ecopipette, CAPP Aero 96 and CAPP Aero 384

Cat. No.	Volume	Use with pipette volumes	Colour
200-FE-100	100 µL	20-200 µL	●
200-FE-200	200 µL	20-200 µL	●
300-FE-300	300 µL	30-300 µL	●
1000-FE-500	500 µL	100-1000 µL	●
1000-FE-1000	1000 µL	100-1000 µL	●
5000-FE-1000	5000 µL	1-5 mL	●
5000-FE-5000	5000 µL	1-5 mL	●
XXX-FE-YYY	Other fixed volume controllers, ecopipette		



C-01



C-04



C-05

Variable Volume Knobs - for CAPP ecopipette, CAPP Aero 96 and CAPP Aero 384

Cat. No.	Volume	Use with pipette volumes	Increment	Colour
2-VZ	0.2-2 µL	0.2-2 µL	0.02 µL	●
10-VZ	0.5-10 µL	0.5-10 µL	0.1 µL	○
20-VZ	2-20 µL	2-20 µL	0.2 µL	●
50-VZ	5-50 µL	5-50 µL	0.5 µL	●
100-VZ	10-100 µL	10-100 µL	1.0 µL	●
200-VZ	20-200 µL	20-200 µL	2.0 µL	●
300-VZ / 1	30-200 µL	30-300 µL	2.0 µL	●
300-VZ	125-300 µL	30-300 µL	2.0 µL	●
1000-VZ-M	100-1000 µL	100-1000 µL	2.0 µL	●
5000-VZ-M	500-5000 µL	500-5000 µL	10.0 µL	●
10000-VZ-M	1000-10000 µL	1000-10000 µL	20.0 µL	●

Pipette Stands

Cat. No.	Description
C-01	Combi stand (2 modules) for up to 3 mechanical pipettes, except CAPPBravo
C-04	Pipette stand for up to 4 CAPP Maestro pipettes or all mechanical pipettes
C-05	Pipette stand for up to 5 CAPP mechanical pipettes, all models
C-06	Pipette stand for up to 6 CAPP mechanical pipettes, except CAPPBravo
00392-1	Block for 4 volume controllers



C-06

Pipette Filters

Cat. No.	Use with pipette	Pipette Volume	Sales Unit
C1000-1-FT	CAPP ecopipette	100-1000 µL	Bag / 50 pcs.
C5000-1-FT	CAPP ecopipette	1-5 mL	Bag / 50 pcs.
B510 mL-25-FT	CAPPBravo	1-5 mL and 1-10 mL	Bag / 25 pcs.



00392-1