

# CAPP Rondo

## GEL CARD CENTRIFUGE

The CAPP Rondo Gel Card Centrifuge is designed for universal compatibility, making it suitable for use with all standard gel cards. It is a very good choice for hospitals, blood banks and laboratories.

### Features

- This excellent centrifuge is capable of holding up to 12 gel cards simultaneously and is compatible with a wide variety of Gel Cards
- Favourite Mode allows users to set up and store up to 99 user-defined programs
- For enhanced safety, the CAPP Rondo Gel Card Centrifuge features a lid lock mechanism that secures the lid during operation
- Safe operation thanks to automatic imbalance detection
- Microprocessor-controlled centrifuge with variable speed and time setting
- The brushless DC motor enables maintenance-free and long-lasting operation
- Compact design with low loading height for easy operation
- Precise RPM setting from 500 to 2200 RPM in steps of 10 RPM
- Safety lid brake, preventing the lid from opening until the rotor has completely stopped
- Functions - last run memory / RPM/RCF speed setting / short spin



### Specifications

<b>Speed</b>	500-2200 RPM / RCF varies	<b>Noise Level</b>	< 65 dB (A)
<b>Speed Accuracy</b>	± 100 RPM	<b>Time Setting</b>	1 to 999 mins and Infinite mode
<b>Motor Type</b>	Brushless DC Motor	<b>Display</b>	Yes
<b>Rotor Capacity</b>	12 Gel cards	<b>Dimension (W x D x H)</b>	381 x 480 x 185 mm
<b>Acceleration Time</b>	35 ± 5 seconds	<b>Weight</b>	8.3kg
<b>Deceleration Time</b>	50 ± 5 seconds	<b>Power Consumption</b>	48 W

## Ordering Information

Cat. No.	Description
CRC-212GC	CAPP Rondo Gel Card Centrifuge 2200 RPM

### Scope of Delivery

1 pc. Gel Card Centrifuge
1 pc. Power Adapter
1 pc. T-Allen Key
Reduction Gel Card Plate
Gel Card Holder
1 pc. User Manual
1 pc. Warranty Card