

Reagent	myDPette Pro
Acetaldehyde	■
Acetic acid 100%	■
Acetic acid 96%	■
Acetic anhydride	■
Acetone	■
Acetonitrile	■
Acetophenone	
Acetylacetone	■
Acetyl chloride	
Acrylic acid	■
Acrylonitrile	■
Adipic acid	■
Allyl alcohol	■
Aluminium chloride	■
Amino acids	■
Ammonia 20%	■
Ammonia 20-30%	■
Ammonium chloride	■
Ammonium fluoride	■
Ammonium sulfate	■
n-Amyl acetate	■
Amyl alcohol	■
Amyl chloride	
Aniline	■
Barium chloride	■
Benzaldehyde	■
Benzol	■
Benzine	
Benzoyl chloride	■
Benzyl alcohol	■
Benzylamine	■
Benzylchloride	■
Boric acid 10%	■
Bromobenzene	■
Bromonaphthalene	■
Butanediol	■
1-Butanol	■
n-Butyl acetate	■
Butyl methyl ether	■
Butylamine	■
Butyric acid	■
Calcium carbonate	■
Calcium chloride	■
Calcium hydroxide	■
Calcium hypochlorite	■
Carbon tetrachloride	
Chloro naphthalene	■
Chloroacetaldehyde 45%	■
Chloroacetic acid	■
Chloroacetone	■
Chlorobenzene	■
Chlorobutane	■
Chloroform	
Chlorosulfonic acid	■
Chromic acid 50%	■
Chromosulfuric acid	■
Copper sulfate	■
Cresol	
Cumene	■

Reagent	myDPette Pro
Cyclohexane	
Cyclohexanone	■
Cyclopentane	
Decane	■
1-Decanol	■
Dibenzyl ether	■
Dichlorobenzene	■
Dichloromethane	
Dichloroacetic acid	■
Dichloroethane	
Dichloroethylene	
Diesel oil	
Diethanolamine	■
Diethyl ether	
Diethylamine	■
1,2 Diethylbenzene	■
Diethylene glycol	■
Dimethyl sulfoxide	■
Dimethylaniline	■
Dimethylformamide	■
1,4 Dioxane	
Diphenyl ether	■
Ethanolamine	■
Ethyl acetate	■
Ethyl alcohol	■
Ethylbenzene	
Ethylene chloride	
Fluoroacetic acid	■
Formaldehyde 40%	■
Formamide	■
Formic acid 100%	■
Glycerol	■
Glycol	■
Glycolic acid 50%	■
Heating oil	
Heptane	
Hexane	
Hexanoic acid	■
Hexanol	■
Hydriodic acid 57%	■
Hydrobromic acid	■
Hydrochloric acid 20%	■
Hydrochloric acid 20-37%	■
Hydrogen peroxide 35%	
Isooctane	
Isoamyl alcohol	■
Isobutanol	■
Isopropanol	■
Isopropyl ether	■
Lactic acid	■
Methoxybenzene	■
Methyl alcohol	■
Methyl benzoate	■
Methyl butyl ether	■
Methyl ethyl ketone	■
Methyl formate	■
Methyl propyl ketone	■
Methylene chloride	
Mineral oil	■

Reagent	myDPette Pro
Monochloroacetic acid	■
Nitric acid 30%	■
Nitric acid 30-70%	■
Nitrobenzene	■
Oleic acid	■
Oxalic acid	■
n-Pentane	
Peracetic acid	■
Perchloric acid	■
Perchloroethylene	
Petroleum	
Petroleum ether	
Phenol	■
Phenylethanol	■
Phenylhydrazine	■
Phosphoric acid 85%	■
Phosphoric acid 85 % +Sulfuric acid 98 %, 1:1	■
Piperidine	■
Potassium chloride	■
Potassium dichromate	■
Potassium hydroxide	■
Potassium permanganate	■
Propionic acid	■
Propylene glycol	■
Pyridine	■
Pyruvic acid	■
Salicylaldehyde	■
Silver acetate	■
Silver nitrate	■
Sodium acetate	■
Sodium chloride	■
Sodium dichromate	■
Sodium fluoride	■
Sodium hydroxide 30%	■
Sodium hypochlorite	■
Sulfuric acid 98%	■
Tartaric acid	■
Tetrachloroethylene	
Tetrahydrofuran	
Tetramethylammonium hydroxide	■
Toluene	
Trichloroacetic acid	■
Trichlorobenzene	
Trichloroethane	
Trichloroethylene	
Trichlorotrifluoro ethane	
Triethanolamine	■
Triethylene glycol	■
Trifluoro ethane	
Trifluoroacetic acid	■
Turpentine	
Urea	■
Xylene	
Zinc chloride 10%	■
Zinc sulfate 10%	■

For your own safety: please observe the general regulations for handling chemicals (e.g. protective clothing and goggles). Read the manual of the dispenser carefully and follow the instructions. Use the dispenser only with regard to its chemical resistance according to this table (status 2017).